

# LED Upgrade



**a|dec**<sup>®</sup>

# Upgrade to Brilliance



**Bright, white, and true.** Replacing your current A-dec light head with the new A-dec LED is quick, easy, and entirely hassle-free. Upgrade to 25% more illuminance while consuming one-fifth of the power needed for halogen bulbs.

**Increase comfort, performance, and reduce your cost of ownership.**

- See colors accurately for better tissue diagnosis
- Gain stunning white light with a color rendering index of 90 at 5,000K
- Reduce shadowing, improve contrast without strain to the eyes
- Optimize ergonomics for the entire dental team
- Position with three axes of rotation and 540° of horizontal motion
- Illuminate without curing with a cure-safe mode at 25,000 lux
- Lower energy costs by 80% and never replace a bulb again
- Treat patients in comfort—without radiating heat or noise

**Am I eligible for an upgrade?** As long as you have a working flexarm on an A-dec light manufactured after 1993, you can upgrade it to an A-dec LED.\*

**Will my new LED be compatible with my A-dec touchpad?** Yes. A-dec touchpads utilizing DCS (Data Communication System) give you the choice to operate the light using the touchpad or local switching. With A-dec DCS, you can power on/off, select from three intensities, activate cure-safe mode, and support auto on/off with chair presets. Without DCS, your new A-dec LED light will have similar control via local switching only.

**How can I learn more?** Visit [a-dec.com/LED](http://a-dec.com/LED) or contact your Authorized A-dec Dealer.

## Ordering Information

Description	Part Number
A-dec LED light head upgrade in White	28.1924.00
A-dec LED light head upgrade in Surf	28.1925.00

\* Compatibility exceptions apply for the *A-dec Performer Light* and some *A-dec 300 lights*. For details, see your Authorized A-dec Dealer or contact an A-dec Customer Service Representative.



**A-dec Headquarters**

2601 Crestview Drive  
Newberg, Oregon 97132 USA  
Tel: 1.800.547.1883 Within USA/CAN  
Tel: 1.503.538.7478 Outside USA/CAN  
Fax: 1.503.538.0276  
[a-dec.com](http://a-dec.com)

*International Distribution Centers*

**A-dec Australia**

Tel: 1.800.225.010 Within AUS  
Tel: +61.(0)2.8332.4000 Outside AUS  
[a-dec.com.au](http://a-dec.com.au)

**A-dec United Kingdom**

Tel: 0800.ADEC.UK (2332.85) Within UK  
Tel: +44.(0).24.7635.0901 Outside UK  
[a-dec.co.uk](http://a-dec.co.uk)