

LUBRICATE & PROTECT

Surgical Milk

Instrument Dip Concentrate

PRE-AUTOCLAVE Instrument lubricant and Protector.

SURGICAL MILK


NEW	Size	Scent
	 4 L	 Unscented

Effective	<i>Bacteriostatic</i>
	<i>Prevents spotting, staining and instrument rust</i>
	<i>Effective in steam, dry or chemical vapour sterilizers</i>
	<i>Routine use reduces cost of instrument repair and replacement</i>
	<i>Will not interfere with the microbial killing efficacy of autoclaving</i>
Convenient	<i>Not sticky, oily or greasy</i>
	<i>1 part solution to 3 parts water can be used for up to 7 days</i>
	<i>Can be used on a variety of metals: steel, carbide, tungsten, stainless steel and aluminum</i>
Safe	<i>Concentrate requires less storage and less waste in our landfills</i>
	<i>Environmentally Friendly. Biodegradable</i>
	<i>Water based and non-toxic</i>

GET THE MOST FROM YOUR INSTRUMENTS



For more information contact your District Manager or call Germiphene 1.800.265.9931. Visit us online at www.Germiphene.com | Email inquiries: info@germiphene.com



got
surgical
milk?

